

Deploy LipidLynxX

Requirement

LipidLynxX is based on FastAPI, thus tool chains using `nginx`, `supervisor`, `gunicorn`, and `uvicorn` is recommended.

- **System installation**

- `sudo apt install nginx supervisor`

- **Python installation**

- Install `anaconda` / `miniconda` and create a virtual enviroment e.g. `envlynx`
 - `conda create -n envlynx python=3.7`
 - `conda activate envlynx`
 - `pip install -r requirement.txt` # do this under LipidLynxX folder
- Install `gunicorn`, and `uvicorn`
 - `pip install gunicorn uvicorn gevent`

Restart linux server before going to following setps.

Setup nginx

- Check if you have `nginx` installed and go to its config folder
 - `cd /etc/nginx/conf.d`
- Create a new config file named `lynx.conf`
 - `sudo nano lynx.conf`
- Edit following content accordingly especially content inside `{}` and save the file.

- ```
server {
 listen 80;
 server_name {example.com};
 access_log /var/log/nginx/example.log;

 location / {
 proxy_pass http://127.0.0.1:8000;
 proxy_set_header Host $host;
 proxy_set_header X-Forwarded-For $proxy_add_x_forwarded_for;
 }
}
```

- `http://127.0.0.1:8000` is default by `gunicorn`

- Restart `nginx`

- `sudo service nginx restart`

## Setup supervisor

---

- Check if you have `supervisor` installed and go to its config folder

- `cd /etc/supervisor/conf.d`

- Create a new config file named `lynx.conf`

- `sudo nano lynx.conf`

- Edit following content accordingly especially content inside `{}` and save the file.

- ```
[program:gunicorn]
directory={path_to_lipid_lynx_folder}
environment=PATH=/home/anaconda3/envs/{envlynx}/bin # or your env path
command=/home/{user_name}/anaconda3/envs/{envlynx}/bin/gunicorn lynx.app:app
-w 4 -k uvicorn.workers.UvicornWorker
autorestart=true
redirect_stderr=true
stdout_logfile={path_to_lipid_lynx_folder}/stdout.log
stderr_logfile={path_to_lipid_lynx_folder}/error.log

[supervisord]
logfile={path_to_lipid_lynx_folder}/supervisord.log
logfile_maxbytes=50MB
```

- if you use miniconda, then replace anaconda3 by miniconda/miniconda3
- `-w 4` means set up 4 workers, you can adjust number of workers based on your server configuration.

- Refresh configs for `supervisor`
 - `sudo supervisorctl reread`
- Restart `supervisor`
 - `sudo service supervisor restart`

Test LipidLynxX

Now you can visit `127.0.0.1` on the PC/server running LipidLynxX and other PCs in the same local network can access LipidLynxX web GUI using the IP of LipidLynxX computer e.g. `192.168.100.101`

Update LipidLynxX

After updating LipidLynxX, you have to restart the web service.

You can reboot the server or run following two commands to restart `supervisor` and `nginx` service.

- `sudo service supervisor restart`
- `sudo service nginx restart`

LipidLynxX should then run with the new version.